




# PAUL<sup>®</sup> GLAUCOMA IMPLANT

PLATE FOR PAUL<sup>®</sup> STABILITY SYSTEM

HOPE FOR  
GLAUCOMA



**AOI**  
Advanced Ophthalmic Innovations



**PAUL® Glaucoma Implant** is a micro-sized tubed aqueous shunt for the PAUL® Stability System which is designed to regulate glaucoma or raised intraocular pressure in the patient's eye.



Drains fluid from the interior of the eye to the exterior of the eye.



Provides a stable & well-formed anterior chamber post-implantation in the immediate post-operation month<sup>1</sup>.




Offers effective & safe glaucoma control for the eye; implanted using the PAUL® Stability System<sup>1</sup>.

<sup>1</sup> Data on file: 22-P2015001-CF02



AOI has received **CE mark** for **Paul® Glaucoma Implant** from the European Commission (EC) and is marketed in United Kingdom, Belgium, France, Germany, Netherlands, Portugal, Spain, Australia, New Zealand, Singapore, Malaysia and Thailand.

## MICRO-SIZED INTERNAL & EXTERNAL DIAMETER TUBE




Smallest tube  
in the market

## PLATE RESERVOIR



Plate surface  
area of 342.1mm<sup>2</sup>


## EFFECTIVE IOP CONTROL



\*Mean IOP of 74 eyes before surgery and at post-operation in 12 months.  
Error bars: 95% Confidence Interval.

\*Data presented in 2020 American Glaucoma Society (AGS).

# INTRODUCING PAUL® STABILITY SYSTEM



- 1 Implant the **PAUL®** Glaucoma Implant
- 2 Apply layer of freeze-dried collagen
- 3 Apply layer of visco-elastic substance



Advanced Ophthalmic Innovations Pte Ltd  
101 Cecil Street #25-04 Tong Eng Building,  
Singapore 069533  
[www.aoi.sg](http://www.aoi.sg) | [quality@aoi.sg](mailto:quality@aoi.sg)

©2018 Advanced Ophthalmic Innovations Pte Ltd.  
All Rights Reserved.  
US Patent No. 9,808,374  
Chinese Patent No. ZL 201580010268.6  
Singapore Patent No. 11201603540P

Advanced Ophthalmic Innovations is the legal manufacturer of  
PAUL® Glaucoma Implant worldwide and sole distributor in Singapore.

APM-P2015001-1805